

# PAXERABURNER

Complete Medical CD/DVD  
Publishing Solution



PaxeraBurner is a complete medical CD/DVD publishing solution for medical imaging departments and scan centers. PaxeraBurner solution enhances both the portability of the imaging studies (no more large bunches of films) and their attainability (viewable on any PC with ordinary specifications).

 **PAXERAMED**  
NEW ERA IN PACS

IHE

CE

FDA





[www.paxeramed.com](http://www.paxeramed.com)

## At a glance

- ✓ Paxera Burner is very easy to use with an intuitive user-friendly touch screen interface and speedy hardware. Paxera Burner also boosts the efficiency of any imaging center by automatic handling of study copying to the optical media. Paxera Burner supports all the CDs and DVDs media available in the market today. This provides a flexibility to handle any type of the studies regardless of the size or the amount of images to copy.
- ✓ Paxera Burner is the one solution for all types of modalities as it supports the DICOM network with the most up-to-date DICOM standards. There is no place that Paxera Burner cannot fit in and do its job flawlessly and efficiently.
- ✓ With Paxera Burner's capabilities of printing customizable labels, you can have different and attractive label designs for your imaging center, aided by the massive customization options for the patient demographics that can be printed on the published discs.

## Features

- ✓ Touch screen for enhanced operation. No need to trouble yourself with keyboards and mice any more.
- ✓ Automatic archiving of studies to CDs or DVDs using state-of-art robot-operated disc publishers.
- ✓ All published CDs/DVDs have built-in DICOM viewer which automatically runs to view the studies when the media is inserted in the drive.
- ✓ Archiving studies to the local storage hard drive.
- ✓ The produced media has PaxeraViewer (DICOM viewer) built-in with auto-run feature.
- ✓ Automatic label printing for the produced CD/DVD with customizable label designs and study information.
- ✓ PACS integration to receive studies from DICOM modalities or PACS workstations.
- ✓ Printing to DICOM or Windows printers with many different printing templates.

## Specifications

- ✓ Publisher model  
Primera BravoPro.
- ✓ Disc Capacity  
Up to 50 discs in standard input/output mode.  
Up to 100 discs using kiosk adapter (included).
- ✓ Recording Speed  
Recording Speeds.  
CD: 40x; DVD: 20x; DVD-DL: 10x.
- ✓ Recordable Formats  
CD: CD-R, CD-RW, CD-Audio (CD-DA), Video-CD, MP3 to CD-Audio, most other industry-standard CD formats.  
DVD: DVD±R, DVD±RW, DVD±DL.
- ✓ Optical Drive  
Latest generation Pioneer drives.
- ✓ Print Method  
Inkjet.
- ✓ Print Resolution  
Up to 4800 dpi; printer driver selectable.
- ✓ Max. Print Width  
4.724" (120mm).
- ✓ Ink Cartridges  
CMY and black monochrome.
- ✓ Colors  
16.7 million and black monochrome.
- ✓ Media Types  
Printable-surface DVD±R and CD-R; standard or water-resistant

